



Fitting Instructions



STEP 1

Wrap Zipper Belt around the chest.

Direction of pull is always across the chest so for right-handed people begin wrapping belt under left arm.

For left-handed people begin wrapping belt under right arm.



STEP 2

Center foam heart over the sternum and fold back excess material and extend loose end of belt to the other side of the belt.



STEP 3

Cut excess material off just in front of the loose plate approximately 1 inch.



STEP 4

After cutting, material should extend past the plate approximately 1 inch.



STEP 5

Direction of pull to achieve chest support is always across the chest.



STEP 6

Once the episode has passed (cough or sneeze), simply push the cord lock button to release support.

* Suspenders are provided to hold zipper belt in place*

zipperbelt.com
(877) 451-3000

