



## Fitting Instructions



Slide PelvicBinder<sup>®</sup> under supine patient. Center the pelvic binder over greater trochanters.



Use gap width between plates as a cutting guide. Cut free end of binder just in front of the plate.



After cutting you should have a 4 to 5 inch gap between the two plates.



After pulling the cord to achieve pelvic reduction the lock will hold the tension at the point you stop pulling. To loosen the binder put slight tension on the cord and press the cord lock button.

## Nursing Considerations

- **Hip circumference range?**  
30" to 61"
- **How long?**  
The PelvicBinders<sup>®</sup> are recommended for temporary pelvic stabilization, 24 to 48 hour use.
- **How often should I check the skin?**  
Daily skin assessments should be performed.
- **Skin assessments?**  
To assess the skin, one health care provider slides his or her fingers under the binder, to hold the pelvis stable, while another health care provider loosens the binder and inspects the skin.

## Intended Use

The PelvicBinder<sup>®</sup> is meant for use by health care professionals to temporarily stabilize fractured bones in the pelvis. The device controls pain and facilitates initial treatment of pelvic trauma.



**pelvicbinder**<sup>®</sup>

**PelvicBinder, Inc.**  
320 CR 4715  
Sulphur Springs, Texas  
75482 USA  
877-451-3000 tel.  
877-451-1820 fax.  
www.pelvicbinder.com



**Emergo Europe**  
Prinsessegracht 20  
2514 AP The Hague  
The Netherlands